# access

## ENDO MOTOR

Scan!

VIDEO

- (en) **Wireless endodontic motor with apex locator**
- (fr) **Moteur d'endodontie sans fil avec localisateur d'apex**
- (it) **Motore endodontico senza fili con localizzatore apicale**
- (es) **Motor de endodoncia inalámbrico con localizador apical**
- (pt) **Motor endodôntico sem fios com localizador apical**
- (de) **Kabelloser Endodontie-Motor mit Apex-Lokalisierer**
- (nl) **Draadloze endodontische motor met apex locator**

---

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## ENDO MOTOR

<b>1. Product description</b>	3
1.1 Components	3
1.2 Scope of application	3
1.3 Contraindications	3
1.4 Main precautions for use	3
1.5 Device safety classification	4
1.6 Main technical parameters and operating environment	4
<b>2. Device installation instructions</b>	4
2.1 Components and accessories	4
2.2 Installation of the contra-angle	5
2.3 Installation of the file	5
2.4 Installation of the apex locator	6
2.5 Installation of the protective sleeves	8
<b>3. Device operating modes and display interface</b>	8
3.1 Motor modes	8
3.2 Apex locator operating mode	9
3.3 Combined motor and apex locator operating mode	9
3.4 Display screens	9
3.5 Terms and definitions	10
<b>4. Instructions for use of the device</b>	11
4.1 Description of the buttons and setting interface	11
4.2 Starting and stopping the motor handpiece	11
4.3 Selecting customized program	11
4.4 Setting the user programme parameters	11
4.5 Selecting file systems	13
4.6 Handpiece functions settings	14
4.7 Protection against torque overload	14
4.8 Limitations of the apex locator	15
<b>5. Troubleshooting</b>	16
<b>6. Cleaning, disinfection, and sterilization</b>	16
<b>7. Storage, maintenance and transportation</b>	18
7.1 Storage	18
7.2 Maintenance	18
7.3 Transportation	19
<b>8. Environmental protection</b>	19
<b>9. After-sales service</b>	19
<b>10. European authorized representative</b>	19
<b>11. Electromagnetic compatibility statement</b>	19

## 1. Product description


The **ENDO MOTOR ACCESS** is suitable for endodontic treatment. It is a cordless endodontic motor with an integrated apex locator. It is used for preparing and enlarging root canals, or as an apex locator device to help determine the root canal working length. It can also be used to prepare canals while monitoring the relative position of the instrument tip within the canal (combined motor and apex locator mode).

This device has the following features:

- Efficient, silent, long-life brushless motor.
- Portable cordless endo motor with combined length determination.
- 360-degree rotation of the contra-angle.
- Real-time feedback technology and dynamic torque control, effectively preventing files from becoming detached.

### 1.1 Components

The device is composed of a charging base, a motor handpiece, a contra angle, a measuring wire, a lip hook, a file clip, a power adapter, a protective silicon cover, etc..

### 1.2 Scope of application

The device must be used in hospitals and clinics by qualified dentists.

### 1.3 Contraindications

- Doctors with pacemakers must not use this device.
- Do not use on patients with pacemakers.
- Do not use on haemophilic patients.
- Use with caution on patients with heart disease, pregnant women, and young children.

### 1.4 ⚠ Main precautions for use

- Please read this instruction manual carefully before using for the first time.
- Do not place this device directly or indirectly near a source of heat. Use and store this device in a reliable environment.
- This device requires special precautions in terms of electromagnetic compatibility (EMC). Do not use this device near fluorescent lamps, radio transmission devices, remote control devices, portable and mobile high-frequency communication devices.
- Do not modify the device. Any modification may violate safety regulations and cause harm to the patient.
- Torque and speed should be set in accordance with the file manufacturer's recommendations.

- Remove the battery from the appliance if it is to be stored for a long time.
- Wireless charging generates heat, so the temperature of the charging base and motor handpiece will rise. It is recommended that the contact time between the handpiece and the base should not exceed 10 minutes during cordless charging.

## 1.5 Device safety classification

- Type of operation mode: Continuous operating device.
- Type of protection against electric shock: Class II equipment with internal power supply.
- Degree of protection against electric shock: B type applied part.
- Degree of protection against harmful ingress of water: Ordinary equipment (IPX0).
- Degree of safety application in the presence of a flammable anaesthetic mixture with air, oxygen, or nitrous oxide: Equipment cannot be used in the presence of a flammable anaesthetic mixture with air, oxygen, or nitrous oxide.
- Applied part: contra angle, lip hook, file clip, touch probe.
- The contact duration of applied part: 1 to 10 minutes.
- Maximum temperature of the applied part : 46.6 °C.

## 1.6 Main technical parameters and operating environment


Power supply mode	Lithium battery, DC 3.7V, 2000 mAh
Power adapter	(Model: ADS-6AM-06N 05050/ UE08WCP-050100SPA) Input: AC 100-240V, 50/60Hz 0,4A Max Output: DC 5V/1A
Torque range	0,4Ncm-5,0Ncm (4mNm ~ 50mNm)
Speed range	100 rpm~1800 rpm
Wireless charging specifications	Frequency range: 112-205KHz Maximum RF output power of the device : 9,46dBuA/m@3m
Environmental parameters	Temperature requirements : 5°C to 40°C Humidity requirements : 30% to 75% Atmospheric pressure requirements : 70 kPa à 106 kPa

## 2. Device installation instructions


### 2.1 Components and accessories

Handpiece

Contra angle


Charging base


Nozzle

Protective silicon cover

Power adapter


## 2.2 Installation of the contra-angle

The contra angle adopts precision gear transmission, and the transmission ratio is 6:1. Use the original contra angle only. Otherwise, the contra angle will be damaged.

### ● 2.2.1 Installation


Align the three pins on the contra-angle with the positioning slots on the motor handpiece. Push the angle-piece horizontally. It will click into place.

The contra-angle is free to rotate, adapting to the root canal in different positions, while monitoring the screen during the operation.

### ● 2.2.2 Removal

Pull the contra-angle out horizontally when the motorised handpiece is not in use.

### ⚠ Precautions for use

- Make sure the contra-angle handpiece is properly connected before starting the motorised handpiece.
- Do not connect or disconnect the contra-angle when the motor is running.


## 2.3 Installation of the file

### ● 2.3.1 Installation

Before starting the device, insert the file into the hole in the head of the contra-angle handpiece. Insert the file while holding down the push button on the contra-angle handpiece.


While pushing, rotate the file back and forth until it is aligned with the locking pins on the contra-angle handpiece.

Once the file is aligned and in position, release the push button to lock the file in the contra-angle handpiece.

### 2.3.2 Removal

Pull on the file while holding down the contra-angle button.


### Precautions for use

- Ensure that the instruments comply with the ISO standard (rotary instrument shanks  $\varnothing 2.334 - 2.350$  mm).
- Connecting and disconnecting files without holding down the push-button can damage the contra-angle handpiece.
- Take care when handling files to avoid injury to your fingers.
- Make sure the file is securely connected and locked before starting the power handpiece.
- Do not disconnect the file while the motor is running.

## 2.4 Installation of the apex locator

a. Push in the button on the file holder, clip the holder onto the metal top of the file, then release the button. In apex locator mode, the file holder must hold the file correctly.


If the connection signal is poor or erroneous, replace the measuring cable.

b. In apex locator mode, it is advisable to install the handpiece in the charging base to obtain a better viewing angle.


- c. The apex locator detects the apical foramen of the canal not the anatomical apex of the tooth. This could explain some of the differences between the apex locator signal and a radiographic image.


Not all conditions are ideal for determining the working length. For the limitations of the apex locator, refer to chapter 4.8.



- d. Connecting the apex locator cables  
Connect the measuring wire to the motor handpiece (USB socket at the rear).

Connect the lip hook to the white socket on the measuring wire.  
Connect the file clip to the black socket on the measuring wire (not necessary in combined motor and apex locator mode).


- e. Test des connexions  
It is strongly recommended to check the quality of the connection before each use.

<p>In apex locator mode, hang the bracket on the lip hook and check that all the indicator bars light up, as shown in the adjacent illustration:</p>	
<p>In combined motor and apex locator mode, touch the lip hook with the file and check that all the indicator bars light up, as shown in the adjacent illustration:</p>	


## 2.5 Installation of the protective sleeves

Before each use and after cleaning and disinfecting the handpiece, fit a disposable insulation sleeve. Insert the insulation sleeve onto the handpiece from the thin end of the handpiece and install it until there are no creases. Once installed, wrap the protective film around the surface of the handpiece. Then clean and disinfect the surface of the handpiece. See section 6 for cleaning and disinfection procedures. After each use, remove insulation sleeve from the thin end of the handpiece.


## 3. Device operating modes and display interface

<p>Please set the operating mode, torque, and speed in accordance with the file manufacturer's recommendations.</p> <p>In conventional motor mode, the torque bar appears on the display (for more information on the torque bar, see section 3.4).</p>	
---	---


### 3.1 Motor modes

<ul style="list-style-type: none"> <li><b>CW Mode (Clockwise continuous rotation mode)</b></li> </ul> <p>The motorised handpiece rotates forward, clockwise rotation.</p>	
<ul style="list-style-type: none"> <li><b>CCW Mode (Anticlockwise continuous rotation mode)</b></li> </ul> <p>The motor handpiece only rotates backwards, anti-clockwise rotation. This mode is used to inject calcium hydroxide and other medicines. In this mode, a continuous double beep sounds continuously.</p>	


<ul style="list-style-type: none"> <li>● <b>ALT Mode</b></li> </ul> <p>F: Front angle, R: Rear angle. The angles of rotation are adjustable.</p>	
<ul style="list-style-type: none"> <li>● <b>ATR Mode (Adaptive torque reversal mode)</b></li> </ul> <p>In this mode, the motor handpiece rotates clockwise and generates a reciprocating movement when the torque on the file exceeds the set limit.</p>	

### 3.2 Apex locator operating mode



<ul style="list-style-type: none"> <li>● <b>EAL Mode (Electronic apex locator)</b></li> </ul> <p>This mode is only for determining the working length of the root canal. In this mode, the motorised handpiece doesn't rotate.</p>	
--	---




### 3.3 Combined motor and apex locator operating mode

When a file is in the canal and the lip hook is in contact with the patient's lip, the unit automatically switches to combined motor and apex locator mode.  
In use, the measuring cable must be connected to the motor handpiece via the USB socket, the white socket is connected to the lip hook and the black socket is inactive.

### 3.4 Display screens

<p><b>Interface de veille</b> <b>Standby interface</b></p> <ol style="list-style-type: none"> <li>Customized program sequence number 0-9</li> <li>Battery consumption</li> <li>Set speed</li> <li>Set torque</li> <li>Operation mode</li> </ol>	
<p><b>Working interface</b></p> <ol style="list-style-type: none"> <li>Set speed</li> <li>Set torque</li> <li>Real time torque</li> <li>Torque display scale</li> </ol>	

<p><b>Canal measurement mode interface (EAL mode)</b></p> <p>a. Apical reference point flash bar b. EAL: Electronic apex locator</p>	
<p><b>Combined motor and apex locator operating interface</b></p> <p>a. File progress indicator bar b. File progress indicator number The numbers 1,0, 2,0, 3,0 (a) and 00-16 (b) do not represent the actual length. They simply indicate the relative position of the file in relation to the apical foramen. c. Apical foramen The number 00 (b) indicates that the file has reached the apical foramen. The numbers -1 and -2 indicate that the file has passed the apical foramen.</p>	
<p><b>Interface for setting the apical reference point</b></p> <p>a. Apical reference point flash bar b. Apical foramen c. Reading of the digital meter</p>	

### 3.5 Terms and definitions


CW	Clockwise rotation, forward direction. Applied to the rotative file.
CCW	Counterclockwise rotation, reverse direction. Applied to the specific files, to inject calcium hydroxide and other solutions.
ALT	The angles of rotation are adjustable. F: Front angle, R: Rear angle.
ATR	Adaptive torque reversal. Up to setting the torque, the motor runs in ATR mode; when the torque is returned to the normal value, the motor turns clockwise.
Forward Angle	Angle of clockwise rotation of the file.
Reverse Angle	Angle of counterclockwise rotation of the file.
EAL	Electronic apex locator. In this mode, the device works like a stand-alone apex.
AP	Apical foramen.
Apical action	The file action when the file tip reaches the flash bar point.
Flash bar position	Indicates the point inside the canal where the apical action is triggered.
Auto start	The file rotation starts automatically when the file is inserted in the canal.
Auto stop	The file rotation stops automatically when the file is taken out of the canal.
Apical slow down	The file slows down automatically as it approaches the apex. Activating in CW and CCW operation mode.
Operation mode	5 operation modes for canal preparing and measuring. CW, CCW, ALT, ATR and EAL.
Speed	Speed of the file rotation.

#### Torque (torque limit/ Trigger)

For CW and CCW modes, the torque value (Torque Limit) triggers the reverse rotation.  
For ATR mode, the torque value (Trigger Torque) triggers the ATR action.

## 4. Instructions for use of the device

### 4.1 Description of the buttons and setting interface

- Preset program selection  
Long press the setting button "P" to access the preset program in standby mode, press the Adjusting button "+"/"-" to select the file system, press the Setting button "P" to access the file number, press Adjusting button "+"/"-" to select file number, then press the Main button to confirm.
- Setting the handpiece functions  
With the handpiece switched off, hold down the «P» adjustment button and press the main button to access the handpiece function setting, press the «P» adjustment button to the target setting, press the '+/-' adjustment button to adjust, then press the main button to confirm.

### 4.2 Starting and stopping the motor handpiece


#### ● 4.2.1 Start

To switch on the device, press the main button. The device will then go into standby mode.  
The screen will display the standby interface of the last programme used. For example:

Standby interface

To switch the motor on from standby mode, press the main button.  
The screen displays the work interface for the last programme used. For example:

Working interface

#### ● 4.2.2 Stop

To stop the motor, press the main button again. The device will return to standby mode. To switch off the device, press the main button while holding down the «P» setting button. In standby mode, the device switches off automatically after 3 minutes. The handpiece switches off automatically when it is placed in the charging base.

### 4.3 Selecting customized program

The motor handpiece has 10 stored programmes (M0-M9) and 5 preset programmes. Press the '+/-' adjustment button to change the sequence number of the custom programme in standby mode. M0-M9 are memory programs for channel preparation and measurement, each memory program has its own parameters such as operating mode, speed and torque, which can be changed.

### 4.4 Setting the user programme parameters

To modify the user programme parameters from standby mode:





1. Press the «P» adjustment button to select the desired parameter.
2. Press the '+/-' adjustment button to change the parameter setting.
3. Press the main button or wait 5 seconds to confirm.




## Precautions for use

- Make sure that the operating mode is suitable before starting the engine.
- All parameters must be set in accordance with the file manufacturer's recommendations.
- Make sure that the settings are correct, otherwise there is a risk of damaging the file.

<p><b>Setting the operating mode</b></p> <p>List of operating modes: CW, CCW, ALT, ATR, EAL (See description of modes in section 3.1)</p> <p>Press the «P» adjustment button in standby mode, then press the '+/-' adjustment button to select the mode.</p> <p>Check that all the parameters in the next level of this operating mode are correct, press the '+/-' adjustment button to select them if they are not.</p>	<p>Operation Mode</p> <p><b>CW</b></p>
<p><b>Setting the working speed</b></p> <p>Press and hold the '+/-' adjustment button to increase or decrease the speed.</p> <p>In continuous rotation modes (CW and CCW), the working speed can be set from 100 rpm to 1800 rpm.</p> <p>In ATR and ALT modes, the working speed can be set from 100 rpm to 500 rpm.</p>	<p>Speed</p> <p><b>250 rpm</b></p>
<p><b>Setting the torque limit</b></p> <p>Press the '+/-' adjustment button to increase or decrease the torque.</p> <p>In CW continuous rotation mode, the torque limit can be adjusted from 0.4 Ncm to 5.0 Ncm depending on the speed set.</p> <p>In ALT mode, the torque limit can be adjusted from 2.0 Ncm to 5.0 Ncm depending on the set speed.</p> <p>In ATR mode, the trigger torque can be adjusted from 0.4 Ncm to 4.0 Ncm.</p> <p>In CCW continuous rotation mode, it is not possible to set a torque limit.</p>	<p>Torque Limit</p> <p><b>2.0 Ncm</b></p>
<p><b>Setting the apical action</b></p> <p>The apical action is triggered when the file reaches the set apical reference point (see Flash bar position). The unit reacts according to the setting you have chosen.</p> <p>Press the '+/-' adjustment button to select.</p> <p>OFF: Deactivates the apical action, the file continues to rotate even if it reaches the reference point.</p> <p>STOP: The unit stops automatically when the file reaches the reference point. It restarts automatically when the file is removed.</p> <p>REVERSE: The device automatically reverses the direction of rotation when the file reaches the reference point. It automatically returns to the original direction of rotation when the file is removed.</p>	<p>Apical Action</p> <p><b>OFF</b></p>
<p><b>Setting the automatic start</b></p> <p>Press the '+/-' adjustment button to select.</p> <p>OFF: Disables automatic start (the main button is required to start the handpiece).</p> <p>ON: The motor starts automatically when the file is inserted into the channel and the file progress indicator displays 2 bars.</p>	<p>Auto Start</p> <p><b>OFF</b></p>
<p><b>Setting the automatic stop</b></p> <p>Press the '+/-' adjustment button to select.</p> <p>OFF: Deactivates the automatic stop (the main button is required to stop the handpiece).</p> <p>ON: The motor stops automatically when the file is removed from the channel.</p>	<p>Auto Stop</p> <p><b>OFF</b></p>

<p><b>Setting the position of the flash bar, apical reference point</b></p> <p>The apex reference point (flash bar) can be adjusted from 2 to AP (apical foramen). The 0.5 scale indicates that the file is very close to the apical foramen. The Apical Action and Apical Slowdown functions are triggered by the apical reference point. Press the '+/-' adjustment button to select the reference point by changing the flash bar.</p>	
<p><b>Apical slow down setting</b></p> <p>ON: the unit slows down to a set final speed as the file tip approaches the apical reference point. OFF: deactivates apical slowing. This function is only available for CW and CCW continuous rotation modes. The final speed must be lower than the nominal speed. Press the '+/-' adjustment button to select.</p>	
<p><b>Forward Angle</b></p> <ul style="list-style-type: none"> <li>In ALT mode, the forward angle can be adjusted from 20° to 400°.</li> <li>In ATR mode, the forward angle can be adjusted from 60° to 400°.</li> </ul> <p><b>Reverse Angle</b></p> <ul style="list-style-type: none"> <li>In ALT mode, the reverse angle can be adjusted from 20° to 400°.</li> <li>In ATR mode, the rear angle can be 20° greater than the front angle.</li> </ul>	  

## 4.5 Selecting file systems

<p>For your convenience, we have preset some file systems. Press the '+/-' button to access the preset program (M0-M9, preset program 1-5), the interface will appear like this.</p>	
<p>Press and hold down the 'P' button to enter the preset program in standby mode, the interface will look like this. Press the '+/-' button to select the file system.</p>	
<p>After selecting the file system, press the 'P' button to enter the file number, press the '+/-' button to select the file number, then press the main button to confirm.</p>	
<p>The preset settings can also be changed from the default settings. If you want to return to the default setting, long press the 'P' button to access the preset program from standby mode, select 'GO-TAPER UNIVERSAL' and press the main button to confirm, the default setting will be restored, switch the handpiece off and on again, the preset program can also restore the default setting.</p>	

## 4.6 Handpiece functions settings

### To set the device parameters:

1. To access the device parameters, with the handpiece motor switched off, hold down the «P» setting button and press the main button.

Software Version


V1.0.0

2. Press the «P» adjustment button to the target setting, then press the '+/-' adjustment button to adjust, finally press the main button to confirm.

<p><b>Auto Power OFF</b> In standby mode, the device switches off automatically once the set timer has elapsed. The timer can be set from 3 minutes to 30 minutes.</p>	<p>Auto Power OFF 5 min</p>
<p><b>Auto Standby Scr</b> The screen returns automatically to the standby interface once the set timer has elapsed. The timer can be set from 3 seconds to 30 seconds.</p>	<p>Auto Standby Scr 30 sec</p>
<p><b>Dominant Hand</b> The device can be set to be left- or right-handed (180° rotation of the display).</p>	<p>Dominant Hand Right</p>
<p><b>Calibration</b> Assurez-vous que le contre-angle d'origine est installé, et Make sure that the original contra-angle is installed and that there is no file before starting the motor calibration. The motor must be calibrated before first use and after lubrication. <b>OFF:</b> No action. <b>ON:</b> Start motor calibration.</p>	<p>Calibration OFF</p>
<p><b>Beeper volume</b> The volume can be set from 0 to 4. Vol. 0: mute.</p>	<p>Beeper Volume Vol.3</p>
<p><b>Restore Defaults (Rétablir les valeurs par défaut)</b> <b>OFF:</b> No action. <b>ON:</b> The device parameters return to the original setting.</p>	<p>Restore Defaults OFF</p>

## 4.7 Protection against torque overload

During operation, if the measured torque load exceeds the torque limit, the motor automatically reverses direction of rotation. The motor returns to the initial operating mode when the torque load falls below the torque limit.

The load value is lower than the preset torque value  
Clockwise rotation

The load value is higher than the preset torque value  
Counterclockwise rotation

Again, the load value is lower than the preset torque value  
Counterclockwise rotation




#### Precautions for use :

- The torque overload protection function is only suitable for CW mode. This function is prohibited in CCW mode, ATR mode.
- When the battery indicator on the motor handpiece shows a low battery capacity, the battery capacity is insufficient for the motor handpiece to reach the torque limit value, i.e. the function will not work properly. Please recharge the battery in good time.
- If the handpiece is constantly charging, the device may stop automatically due to the overheating protection. In this case, switch off the motor handpiece for a while until the temperature drops.

### 4.8 Limitations of the apex locator

Not all conditions are ideal for assessing root canal length. An accurate signal can only be obtained if the root canal meets the conditions listed below.

<p><b>Root canal with a wide apical foramen</b> A root canal with an exceptionally wide foramen due to injury or incomplete development can interfere with the electrical signal. The results may indicate a shorter length than the actual one.</p>	
<p><b>Root canal with liquid flowing from the canal entrance</b> A root canal in which blood or other liquid flows from the canal entrance and encounters the gum tissue can disrupt the electrical signal, making it impossible to obtain an accurate measurement. Wait until the flow has stopped completely, clean the canal, and then proceed with the measurement.</p>	
<p><b>Broken crown</b> If the crown is broken and part of the gum tissue penetrates the cavity surrounding the canal entrance, contact between the gum tissue and the file may disrupt the electrical signal. In this case, cover the tooth with a suitable material to isolate the gum tissue.</p>	
<p><b>Fractured tooth</b> A fractured tooth can disrupt the electrical signal.</p>	

<p><b>Root canal previously filled with gutta-percha</b> Gutta-percha debris can interfere with the electrical signal. The gutta-percha must be completely removed to eliminate the insulating effect.</p>	
<p><b>Crown or metal prosthesis touching gum tissue</b> Contact between the prosthesis and the file may disrupt the electrical signal. In this case, widen the opening at the top of the crown so that the file does not touch the metal prosthesis before taking a measurement.</p>	
<p><b>Extremely dry root canal</b> A dry root canal can interfere with the electrical signal. In this case, try moistening the canal with saline solution.</p>	

## 5. Troubleshooting

Fault	Possible cause	Solution
The motor handpiece does not rotate.	It is in EAL mode, as it is only for canal measurement.	Change to CW, CCW, ALT or ATR mode.
There is a continuous beep when the motor handpiece is running.	The continuous beep sound is indicating that the motor handpiece is under CCW mode.	Stop the motor handpiece and change the operation mode to CW Mode.
The contra angle calibration failed.	The calibration failure can be caused by a strong resistance of the contra angle.	Clean the contra angle and recalibrate it after the oil injection.
The endurance time becomes shorter after charging.	The battery capacity is decreasing.	Please contact the local distributor or manufacturer.
No sound.	Beeper Volume set to 0. Vol.0: Mute.	Set the Beeper Volume to 1,2,3.
The rotating file is stuck in the root canal.	Incorrect specification setting. File load torque set too high.	Choose the CCW Mode, start the motor handpiece, and take the file out.

## 6. Cleaning, disinfection, and sterilization

For reasons of hygiene and health safety, the motorised handpiece, AC adapter and charging base must be cleaned and disinfected before each use to avoid contamination. In addition, the contra-angle handpiece, lip hook, file clip, protective silicone cover and probe must be cleaned, disinfected, and sterilised before each use to prevent contamination. This applies to the first use, as well as all subsequent uses.

### General recommendations

- Before and after each use, all objects that have been in contact with infectious agents must be cleaned using towels impregnated with a disinfectant.
- Only use a disinfectant solution that has been approved for effectiveness (VAH/DGHM list, CE mark, FDA and Health Canada approval) and in accordance with the disinfectant solution manufacturer's DFU.
- Do not immerse the contra-angle handpiece in a disinfectant solution or in an ultrasound bath.
- Do not use bleach, detergents or chloride-based disinfectants.

- For your own safety, please wear personal protective equipment (gloves, goggles, mask).
- It is the user's responsibility to ensure that the equipment is sterile.
- Water quality must comply with local regulations, particularly for the final rinsing stage or with a washer-disinfector.
- Refer to the manufacturer's instructions for use for sterilising endodontic files.
- The contra-angle handpiece must be lubricated after cleaning and disinfection, but before sterilisation.
- Do not sterilise the motorised handpiece, mains adaptor and base.

## 6.1 Cleaning, disinfection, and sterilization of the handpiece, the adapter and base

- The handpiece, adapter and base cannot be cleaned and disinfected using automatic equipment. Manual cleaning and disinfection are necessary.
- Cleaning and disinfection must be carried out within 10 minutes of use.

Step	Process
1. Preparation	Retirez la pièce à main, l'adaptateur et la base de la zone de travail. Retirez les accessoires (contre-angle, crochet à lèvres, porte-lime, sonde de contact) de la pièce à main et de la base.
2. Manual cleaning	Wet a soft, clean cloth with distilled or demineralised water at room temperature and wipe the surface of the components thoroughly 3 times. Replace the cloth with a dry, soft, lint-free cloth after each wipe. If any visible contaminants remain, wipe until they have been removed.
3. Manual disinfection	Soak a soft, dry cloth in 75% alcohol and wipe all surfaces of the handpiece, adapter and base for at least 3 minutes. Then wipe them dry with a soft, lint-free cloth.
4. Control	Inspect the motorised handpiece, mains adaptor and base and remove any defective components. If any component is soiled, repeat the cleaning and disinfection process.
5. Storage	Keep the handpiece, charger, base and other components in a clean place.

## 6.2 Cleaning, disinfection, and sterilization of the accessories (contra-angle, lip hook, file clip, touch probe)

- The use of strong detergents and disinfectants (alkaline pH>9 or acid pH <5) will reduce the life of the products. In this case, the manufacturer declines all responsibility.
- The accessories have been designed for many sterilisation cycles. However, with each new procedure to prepare for use, thermal and chemical stresses cause the products to age (maximum number of sterilisations = 250 times).
- Cleaning, disinfection and sterilisation must be carried out within 30 minutes of use.

Step	Process
1. Preparation	Products must be processed in a disassembled state. Remove the accessories and disconnect the file and contra-angle from the handpiece and base. Carry out manual pre-cleaning until the handpiece is visually clean, to remove most of the dirt. The water temperature should not exceed 40°C.
2. Automated cleaning / disinfection / drying in a washer-disinfector	The washer-disinfector must meet the requirements of standard ISO 15883. Carefully place the accessories in the washer-disinfector. Ensure that the accessories cannot move or encounter each other. The cycle: Ao value >3000 or, at least 5 min at 90°C/194°F The solution used may be pure water, distilled water, deionised water, or a multi-enzyme solution.
3. Drying	Can be washer disinfected. If necessary, additional manual drying can be carried out using a lint-free cloth. Blow sterile compressed air into the instrument cavities (maximum pressure 3 bar). Accessories can also be dried in a medical drying cabinet or oven (80°C for 15 minutes).
4. Control	Inspect the accessories and sort out those with defects. Accessories considered soiled must be cleaned and disinfected again. If the accessories are damaged, replace them. Lubricate the contra-angle with a suitable spray before packing Refer to chapter 7.2.
5. Packaging	Accessories must be quickly packed in a sterilisation bag or sterile box). Use a sterilisation bag which complies with standard EN ISO 11607-1 and is resistant to temperatures of up to 137°C. Avoid contact with dissimilar metals during packaging.

6. Sterilization	Sterilization of accessories using a fractionated pre-vacuum steam sterilisation process (in accordance with EN 285/EN 13060/EN ISO 17665), considering the requirements of the respective countries. Minimum requirements: at least 4 minutes at 132°C/134°C (in the EU: 5 minutes at 134°C, in the USA: 4 minutes at 132°C).
7. Storage	Store the accessories in the sterilisation packaging in a dry, clean environment. Sterility cannot be guaranteed if the packaging is opened, damaged or damp. Check the packaging and the corner before use (integrity of the packaging, absence of moisture and shelf life). Shelf life should not exceed 7 days.

## 7. Storage, maintenance, and transportation

### 7.1 Storage

- The equipment must not be stored with toxic, corrosive, flammable or explosive items.
- This equipment should be stored in a room where the relative humidity is between 10% and 93%, the atmospheric pressure between 70kPa and 106kPa, and the temperature between -20°C and +55°C.

### 7.2 Maintenance


#### • 7.2.1 Calibration

Calibrate after replacing or lubricating the contra-angle (see section 4.6).

#### • 7.2.2 Lubrication of the contra-angle

The contra angle needs to be lubricated after cleaning and disinfection, but before sterilization.

1. Screw the oil injection nozzle onto the oil bottle (approximately 1 to 3 turns).
2. Insert the nozzle into the end of the contra-angle.
3. Fill the angle-piece with oil until the oil flows from the head of the angle-piece.
4. Hold the angle-piece firmly so that it is not propelled by the pressure.
5. Position the angle-piece vertically for at least 30 minutes to allow the excess oil to drain by gravity through the tip.

#### Precautions

- Only the original oil injection nozzle may be used to lubricate the contra-angle.
- The motor handpiece must not be filled with oil.

#### • 7.2.3 Battery charging

The motorised handpiece has a built-in rechargeable lithium battery.

1. Insert the mains adapter plug into the power socket on the charging base and check that they are correctly connected.
2. Then insert the motor handpiece into the charging base (the motor handpiece must be correctly aligned with the charging base).  
When the blue light on the charging base flashes, charging is in progress.  
When the handpiece is fully charged, the blue light on the charging base stays on.
3. After charging, unplug the adapter.

#### • 7.2.4 Battery replacement

1. Replace the battery if it seems to run down sooner than expected.
2. Switch off the device.
3. Using tweezers or a screwdriver, open the rubber cover and remove the screw.
4. Remove the battery cover.
5. Remove the battery to be replaced and disconnect the connector.
6. Connect the new battery and place it in the motor handpiece.
7. Replace the cover and its screw.

## 7.3 Transportation

- Prevent excessive impact and vibration during transportation.
- During transportation, the product should not be mixed with dangerous goods.
- Avoid exposure to the sun or rain or snow during transportation.

## 8. Environmental protection

Please dispose according to the local policies.

## 9. After-sales service

As of the date of sale, if the equipment does not work normally due to quality problems, our company shall be responsible for the maintenance with the use of the warranty card. The warranty period and warranty scope shall refer to the product warranty card. This product does not contain parts suitable for self-repair, any maintenance should be carried out by designated professionals or a special maintenance shop. If you need to repair the equipment components, the circuit diagram, component list, correction details and other information can be obtained by contacting the manufacturer.

## 10. European authorized representative

**MedNet EC-REP C IIb GmbH**  
Borkstrasse 10 · 48163  
Muenster · Germany

## 11. Electromagnetic compatibility statement

The device has been tested and approved in accordance with EN 60601-1-2 for EMC. This does not guarantee that this device will not be affected by electromagnetic interference. Avoid using the device in a high electromagnetic environment.

### 11.1 Technical description concerning electromagnetic emission

Table 1: Guidance and declaration of the manufacturer - electromagnetic emissions

Guidance and declaration of the manufacturer - electromagnetic emissions		
The Endo Motor ACCESS is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - Guidance
RF emissions CISPR 11	Group 1	The Endo Motor ACCESS The model Ai-Motor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Endo Motor ACCESS is suitable for used in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

### 11.2 Technical description concerning electromagnetic immunity

Table 2 : Guidance and declaration of the manufacturer - electromagnetic immunity

The Endo Motor ACCESS is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance

Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ±4, ±8, ±15kV air	± 8 kV contact ± 2 kV, ±4, ±8, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transients/bursts IEC 61000-4-4	± 2kV for power supply lines ±1kV for input/output lines	± 2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Overvoltage IEC 61000-4-5	±0,5, ±1kV line to line ±0,5, ±1, ±2kV line to earth	±0,5, ± 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations in the power supply IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle <5 % UT (>95 % dip in UT) for 1 cycle 70% UT (30 % dip in UT) for 25 cycles <5% UT (>95 % dip in UT) for 250 cycles	<5 % UT (>95 % dip in UT) for 0,5 cycle <5 % UT (>95 % dip in UT) for 1 cycle 70% UT (30 % dip in UT) for 25 cycles <5% UT (>95 % dip in UT) for 250 cycles	The quality of the mains supply should be that of a typical commercial or hospital environment. If the user requires continuous operation during power cuts, it is recommended that the device is powered by an inverter or battery.
Power frequency magnetic field (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	Magnetic field levels from the mains voltage should be those of a typical commercial or hospital environment.

Note: UT is the AC mains voltage before the test level is applied.

**Table 3 : Guidance and declaration - electromagnetic immunity with regard to radiofrequencies emitted and radiated**

Guidance and declaration of the manufacturer - electromagnetic immunity			
The Endo Motor ACCESS is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V (valid value) 150 kHz to 80MHz	3 V (valid value)	Portable and mobile RF communications equipment should be used no closer to any part of the Endo Motor ACCESS including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80MHz to 2,7GHz	3 V/m	Recommended separation distance $d=1,2 \times P^{1/2}$ $d=2 \times P^{1/2}$ $d=1,2 \times P^{1/2}$ 80 MHz to 800 MHz $d=2,3 \times P^{1/2}$ 800 MHz to 2,7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range, b) interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz end 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Endo Motor ACCESS is used exceeds the applicable RF compliance level above, the Endo Motor ACCESS should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Endo Motor ACCESS. Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

**Table 4 : Recommended separation distances between portable and mobile RF communications equipment and the Endo Motor ACCESS**

The Endo Motor ACCESS is intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{p}$	80 MHz to 800 MHz $d=1.2\sqrt{p}$	800 MHz to 2,5 GHz $d=2.3\sqrt{p}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

**NOTE 1** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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