# access

## EASY PACK

- ⓔn Obturation system for root canal filling
- ⓕr Système d'obturation du canal radicaire
- ⓔs Sistema de obturación del conducto radicular
- ⓖt Sistema de obturação de canais radiculares
- ⓔt Sistema di otturazione del canale radicolare
- ⓔe Obturationssystem für Wurzelkanalfüllungen
- ⓔl Obturatiesysteem voor wortelkanaalvulling

---

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02 - 11  
12 - 21  
22 - 31  
32 - 41  
42 - 51  
52 - 61  
62 - 71

- 1. Product introduction**
  - 1.1 Package contents
  - 1.2 Scope of application
  - 1.3 Contraindications
  - 1.4 Equipment safety classification
  - 1.5 Main technical parameters and operating environment
- 2. Product installation and function description**
  - 2.1 Schematic diagram of Central Unit
  - 2.2 Schematic diagram of components
  - 2.3 Power adapter installation
  - 2.4 Use of each function button
  - 2.5 Functional operation description
  - 2.6 Assembly and disassembly of working tips
- 3. Product use instructions**
- 4. Precautions**
- 5. Troubleshooting**
- 6. Cleaning, disinfection, and sterilization**
- 7. Storage, maintenance, and transportation**
  - 7.1 Storage and maintenance
  - 7.2 Transportation
- 8. Environmental protection**
- 9. After-sales service**
- 10. Electromagnetic compatibility statement**
- 11. Symbol description**

## 1. Product introduction

The **Easy Pack** obturation system is suitable for root canal filling during root canal treatment. It is used to soften the gutta percha by heating it, and then fill it into the prepared root canal.

The equipment features:

- Wireless design, small and easy to hold;
- Clear, simple and easy-to-operate LED screen;
- Precise temperature control, five preset temperatures can be selected: 90 °C, 150 °C, 180 °C, 200 °C, 230 °C;
- Safety protection mechanism, stopping heating after 5S automatically.

### 1.1 Package contents

Name	Quantity
1. Central unit	1
2. Charging base	1
3. Power adapter	1
4. Type-C charging wire	1
5. Working tip	4
6. Instructions	1

### 1.2 Scope of application

The **Easy Pack** obturation system is suitable for filling the root canal during root canal treatment.

### 1.3 Contraindications

- Prohibited for patients with known allergy to natural latex, stainless steel, silver, copper and other metals.
- Prohibited for patients with hemophilia.
- Prohibited for patients with pacemakers.
- Prohibited for use by physicians with pacemakers.
- Use with caution in patients with heart disease, pregnant women and young children.

### 1.4 Equipment safety classification


- Classification by operation mode: short-time operation.
- Classification by electric shock protection type: class II equipment.
- Classification by electric shock protection degree: type B applied part.
- Protection against ingress of fluids: ordinary equipment (IPX0).
- Not to be used with a flammable anesthetic gas mixed with air or a flammable anesthetic gas mixed with oxygen or nitrous oxide.

### 1.5 Main technical parameters and operating environment

Size	Central unit: 160 mm x 24 mm x 24 mm
Weight	About 73 g
Power supply mode	Lithium battery, DC3.7V+ 10%, 1900 mAh
Power adapter	Input: 100-240V, 50 / 60Hz, 0.4 A Output: DC5V/1A
Fluid penetration protection	IPX0
Classification by electric shock protection type	Class II, internal power supply equipment
Classification by electric shock protection degree	Class B
Optional temperature	90°C-150°C-180°C-200°C-230°C
Power dissipation	<3W
Service environment	Temperature requirements: + 5°C ~ + 40°C
Storage / transport environment:	Temperature requirements: -20°C ~ + 55°C Humidity requirements: 10% -93% Atmospheric pressure: 500 hPa~1060 hPa

## 2. Product installation and function description

### 2.1 Schematic diagram of Central Unit

### 2.2 Schematic diagram of components


1. Power button
2. LED display screen
3. Heating button
4. Central unit

5. Type-C charging cable
6. Charging base
7. Power adapter
8. Working tip

### 2.3 Power adapter installation

This product uses contact charging. When the power is low, please put the unit on the charging base for charging.

Connect the power adapter and base to both ends of the type-C charging cable as shown in Figures (a) and (b). When the central unit is correctly connected to the charging base, the charging base display light changes from blue to flashing green; after full charge, the light turns to steady green.

## 2.4 Use of each function button

### ● 2.4.1 Power button

When the device is off, long press the power button to switch it on.

When the device is on, long press the power button to switch it off.

When the device is on, short press the power button to set the temperature, as shown in the figure below.

### ● 2.4.2 Heating button

When the device is on, long press the heating button to raise the temperature.

## 2.5 Functional operation description

### ● Main interface

1. After powering on, the main interface displays the default temperature value (150 °C), STANDBY state (STANDBY) and real-time battery power as shown in the figure below.

2. Battery icon: As shown in the figure below, the battery power display icon is dynamic during battery charge.


				
80-100% electric charge	60-80% electric charge	40-60% electric charge	20-40% electric charge	0-20% charge, the battery is depleted. Please charge the device immediately.

**Note :** This icon indicates the remaining battery power. During gutta percha filling, the battery power will decrease.

### ● 2.5.2 Heating interface

Long press the heating button for heating. The display screen displays the real-time heating temperature, the heating state (HEATING) and the remaining battery power.


### ● 2.5.3 Cooling interface

When the heating button is released or the heating time is more than 5S, the system starts cooling automatically, the display screen displays the real-time cooling temperature and cooling state (COOLING).

## 2.6 Assembly and disassembly of working tips

### ● 2.6.1 Assembly of the working tip

When assembling the working tip, as shown in Figure (a), insert the hexagonal end of the working tip (Figure c) into the hexagonal socket of the central unit (Figure b).

- **2.6.2 Disassembly of the working tip**

Before removing the working tip, make sure that it has cooled down, then pull it off.

### 3. Product use instructions

Working tip selection

The working tip should be selected according to the actual situation of the patient. Each tip size and taper is marked with a different color ring. Details are shown in the following table:

Working tip size:

Type	Specification	Taper
SN-F04	F	04
SN-FM06	FM	06
SN-M08	M	08
SN-ML10	ML	10

### 4. Precautions

- As it is the case with all electronic devices, this instrument causes electromagnetic interference and should not be used in patients with pacemakers.
- During operation and cleaning, pay attention to the direction and position of the gutta percha injection needle to avoid injuring patients and operators.
- Pay full attention to whether the gutta percha injection needle is loose, hot or makes noise. If you detect any abnormal phenomenon, please stop using the device and contact the local dealer or manufacturer.
- Keep the device safe from impacts and especially from dropping.
- When removing the gutta percha injection needle and heat insulation protection cover, turn off the power first in order to avoid any accidental starting and injury that might result from touching the trigger accidentally.
- When the battery power icon shows low power, please charge in time using the original power adapter.
- Please use the original accessories, those from other manufacturers may not work with this equipment or cause equipment damage.
- Please do not disassemble or repair the equipment without permission, otherwise the warranty will no longer be valid.
- This equipment should not be placed in wet places or be in contact with any liquid.
- Do not expose the equipment to a direct or indirect heat source. The equipment should be operated and kept in a safe environment.

### 5. Troubleshooting

Fault	Possible cause	Solution
The device cannot power on normally.	Insufficient battery power.  Short circuit of the charging interface.  The power button pressing time is too short.	Charge the device connecting it to a power supply.  Remove the foreign bodies at the interface, wipe it clean and reconnect it to the power supply.  Long press the power button to start the device.
The charging base display light is not on.	Power adapter specification mismatch.	Reconnect after checking.
With the power adapter plugged in, the device is not charged.	Not correctly plugged in.  Power supply is damaged or the specifications do not match.  The charging base contacts are dirty.	Pull out and reconnect.  Contact your dealer.  Wipe the charging base contacts with alcohol.
With the battery fully charged, the service cycle becomes shorter.	Reduced battery capacity	Return the device to the manufacturer to replace the battery.

The display screen displays TIP ERROR when heating	The working tip is not inserted correctly.	1. Insert the working tip.
	The working tip is damaged.	2. Replace it with a new tip. 3. Contact the after-sales service.

## 6. Cleaning, disinfection and sterilization

### 6.1 Working tip cleaning and disinfection

Step	Process
1. Rinse	Wash the gutta percha needle, silver needle bending device, heat insulation protection cover with running water for 2 min to remove the surface contaminants.
2. Wipe	Wet a soft clean cloth with a cleaning agent and thoroughly wipe the surface of the product 5 times. Replace the soft cloth after each wipe. If there are still visible pollutants remaining, wipe repeatedly until there are no visible pollutants.
3. Brush wash	Brush the product thoroughly with an instrument brush for 3 min.
4. Soak	Immerse the product in the cleaning agent for 5 min.
5. Rinse	Wash the product with purified water for 2 min to remove any residual detergent from its surface.
6. Dry	Use a dry soft absorbent cloth to eliminate any residual water trace on the product surface.

### 6.2 High-temperature and high-pressure sterilization

Sterilization	After cleaning, put the working tip in a disposable sterilization bag and carry out sterilization at a temperature of 134°C during 18 minutes.
---------------	--

## 7. Storage, maintenance and transportation

### 7.1 Storage and maintenance

This equipment should be kept away from any source of heat and installed or kept in a cool, dry and ventilated place.

The equipment should not be stored with toxic, corrosive, flammable and explosive articles.

The equipment should be stored in an environment with a relative humidity of 10% -93%, an atmospheric pressure of 500 hPa ~ 1060 hPa and a temperature of -20°C ~ + 55°C.

When the equipment is not in use, turn off the power switch and unplug the charging base.

Using the equipment will shorten the battery life. Please charge in time when the power is low.

When the equipment is not in use for a long time, it should be charged once a month for 1 hour.

### 7.2 Transport

Prevent excessive impact and vibration during transportation.

During transportation, the product should not be mixed with dangerous goods.

Avoid exposure to the sun or rain or snow during transportation.

## 8. Environmental protection

The equipment does not contain any harmful ingredients and can be treated or disposed of in accordance with relevant local policies.

## 9. After-sales service

As of the date of sale, if the equipment does not work normally due to quality problems, our company shall be responsible for the maintenance with the use of the warranty card. The warranty period and warranty scope shall refer to the product warranty card. This product does not contain parts suitable for self-repair, any maintenance should be carried out by designated professionals or a special maintenance shop. If you

need to repair the equipment components, the circuit diagram, component list, correction details and other information can be obtained by contacting the manufacturer.

## 10. Electromagnetic compatibility statement

The **Easy Pack** downpack obturation device meets the relevant requirements of the YY9706.102-2021 standard on EMC, and must be installed and used according to the EMC information specified in this manual. Portable and mobile RF communication equipment may affect this device.

Except for cables (transducers) sold as spare parts for internal components, the use of accessories and cables (transducers) other than those specified for the device may result in increased emissions or decreased immunity of the equipment or system.

The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.

### 10.1 Guidance and manufacturer's declaration - electromagnetic emissions

Guidance and manufacturer's declaration - electromagnetic emissions		
The Easy Pack Downpack Obturation Device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment		
Emissions test	Compliance	Electromagnetic environment - Guidance
Radio-frequency emissions GB-4824	Group 1	The Easy Pack Downpack Obturation Device uses RF energy only for its internal function. Therefore, its RF emissions are low and are not likely to cause any interference in nearby electronic equipment.
Émissions de radiofréquences GB-4824	Class B	The Easy Pack Downpack Obturation Device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes.
Radio-frequency emissions GB-4824	Class A	
Harmonic emissions GB -17625.1	Complies	

### 10.2 Guidance and manufacturer's declaration - electromagnetic immunity

Guidance and manufacturer's declaration - electromagnetic immunity			
The Easy Pack Downpack Obturation Device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.			
Immunity test	IEC-60601 test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge GB/T 17626.2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst GB/T 17626.4	± 2kV for power supply lines	± 2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge GB/T 17626.5	± 1 kV line to line ± 2 kV line to earth	± 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines GB/T 17626.11	<5 % UT (> 95% dip in UT) for 0.5 cycle  40 % UT (60% dip in UT) for 5 cycles  70 % UT (30% dip in UT) for 25 cycles  <5 % UT (> 95% dip in UT) for 5 sec.	<5 % UT (> 95% dip in UT) for 0.5 cycle  40 % UT (60% dip in UT) for 5 cycles  70 % UT (30% dip in UT) for 25 cycles  <5 % UT (> 95% dip in UT) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Easy Fill obturation system requires continued operation during power mains interruption, it is recommended that the Easy Fill obturation system be powered from an uninterruptible power supply or battery.
Power frequency (50/60Hz) magnetic field GB/T 17626.8	3 A/m	3 A/m (50/60Hz)	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the AC mains voltage prior to application of the test level.

## 10.3 Guidance and manufacturer's declaration - electromagnetic immunity

### Guidance and manufacturer's declaration - electromagnetic immunity

The Easy Pack Downpack Obturation Device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
Conducted RF GB/T 17626.6	3 V (valid value) 150 kHz-80MHz	3 V (valid value)	Portable and mobile RF communications equipment should be used no closer to any part of the Easy Fill obturation system including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance $d = 1.2$ $d = 1.2 \sqrt{P}$ 80MHz~800MHz $d = 2.3 \sqrt{P}$ 800MHz~2.5GHz  where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range, b) interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF GB/T 17626.3	3 V/m 80MHz-2.5GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Easy Fill obturation system is used exceeds the applicable RF compliance level above, the Easy Fill obturation system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Easy Fill obturation system. Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

## 10.4 Recommended separation distances between portable and mobile RF communications equipment and Easy Fill

The Easy Pack Downpack Obturation Device is intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter in meter		
	150 kHz~80MHz $d=1.2\sqrt{p}$	80MHz~800MHz $d=1.2\sqrt{p}$	800MHz~2.5GHz $d=2.3\sqrt{p}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

## 11. Symbol description

	Refer to instructions		Warning: Read the instructions before use
	Type B applied part		Medical device
	The product complies with the WEEE directive, and the equipment must be treated as municipal solid waste when disposed of		Product serial number
	Manufacturer		Date of manufacture
	Indoor use		Upward
	Keep dry		Fragile, handle with care
	UDI		

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EQUIPMENT CO., LTD**