

access

NEW

ENDO MOTOR

Wireless endodontic motor with integrated apex locator

- **Facilitates root canal preparation**
- Real-time feedback technology and dynamic torque control
- Automatic stop, slow down or reverse rotation when the file approaches or reaches the apical foramen
- Equipped with **27** pre-registered systems and **10** customisable programmes

Autostart mode

The motor starts up when it enters into the canal

ATR mode

When the torque limit is reached the motor changes direction of rotation



ENDO MOTOR

Brushless motor !

- **360°** rotation of the contra-angle
- Precise continuous rotation and alternate movement
- Clear, easy-to-use **LED screen** and convenient **wireless** charging!
- **Speed:** 100 – 1800 rpm
- Torque range: 0.4 – 5 Ncm
- Weight of the handpiece with the contra-angle: **158 g**
- Frequency range: 112 – 205 KHz



Motor
2-year warranty

Contra-angle
1 year warranty

Supplied with

Instructions for use, quick guide, charging base, handpiece, contra-angle, measuring cable, lip hook, file holder, mains adaptor, protective silicone cap, spray nozzle, touch probe.

ENDO MOTOR

3166410

