

# access

## EASY FILL

Gutta-percha injector designed for obturation of the 2/3 coronal of the canal



- Facilitates injection of gutta-percha
- 4 preset temperatures: 100°C, 150°C, 180°C and 200°C
- Clear and easy-to-use LED screen
- Cordless gun-shaped obturation system!

Supplied with

6 EASY FILL TIPS (2 x 25G-24MM, 2 x 25 G-28MM, 2 X 24G-24MM), 1 tip bending tool and 3 thermal insulation sleeves.

EASY FILL

3153026

## EASY FILL

### TIPS

x6

134°C



Sterilisable at 134°C for 18 minutes



#### EASY FILL TIPS

25G (0,45MM) 24MM	3153027
25G (0,45MM) 28MM	3153028
24G (0,55MM) 24MM	3153029
23G (0,65MM) 24MM	3153030
23G (0,65MM) 28MM	3153031

# access

## EASY PACK

Hot condenser  
for root canal filling

- Used to soften gutta-percha then compact it
- 5 preset temperatures: 90°C, 150°C, 180°C, 200°C and 230°C
- Cordless small obturation system !



Supplied with  
4 assorted EASY PACK TIPS.

EASY PACK

3153032

## EASY PACK TIPS

x4



### EASY PACK TIPS

04/40	3153033
06/55	3153034
08/55	3153035
10/55	3153036

134°C



Sterilisable at 134°C  
for 18 minutes